


AHCCCS Health Information Exchange and Electronic Health Record Utility Phase I Technical Overview The HIEHR Utility Project April 2008




Conceptual View of AHCCCS HIE















HIE Architecture






Viewer Login






Patient Search: By Name






AHCCCS Test System Results Displayed





AHCCCS Test System Patient Search: By AHCCCS ID






Results Displayed






AHCCCS Test System Patient Details Displayed: Filter Records






AHCCCS Test System Patient Details Displayed: Filter Records





AHCCCS Test System Select Multiple Patients






AHCCCS Test System Patient Details Displayed

		212200 CONT				
Patient Records - Microsoft Internet Explorer						
Elle Edit <u>V</u> iew F <u>a</u> vorites <u>T</u> ools <u>H</u> elp					Links »	
AL/GGGS						
<u>Viewer Login</u> > <u>Patient Search</u> > Patient Records Log					ogged In: Kelly, Mary <u>Loq</u>	
					N N	3
Filter Criterion						<u>elp</u>
Select your Record Ty Select Dates: O Tim	1.55 (4.55)					
O Date Range From: (mm/dd/yyyy) IIII To: (mm/dd/yyyy) IIII						
	78					
					Apply Filter Clear	
PATIENT RECORDS						
						=
Last Name: SMITH First Name: BOB MI; Date of Birth: Gender: M AHCCCS Id: A95195195 Last 4 digits of SSN: Address: 123 KELLER LANE PHOENIX AZ 86413					05/15/1970 0444	
Sel Record Type	Date of Service	Source	Category	Record Descriptor		
Discharge Summary	05/11/2002	AHCCCS	Α	Discharge Summary		Details for both
Lab Report	11/01/2001	AHCCCS	А	Lab Report		Patients displayed
Medication History	05/11/2002	AHCCCS	A	Medication History		
Last Name: SMITHE First Name: ROBERT MI: Date of Birth: 05/15/1970 Gender: M AHCCCS Id: A95195196 Last 4 digits of SSN: 0444 Address: 123 KELLER LANE PHOENIX AZ 86413						
Sel Record Type	Date of Service	<u>Source</u>	Category	Record Descriptor	1.0	-
Discharge Summary	07/21/1992	AHCCCS	Α	Discharge Summary		
Lab Report	09/27/2007	AHCCCS	A	Lab Report		
Medication History	03/01/2008	AHCCCS	A	Medication History		<u>♥</u>
Pone					S Local intranet	




Patient Details Displayed: Filter On 🛂 Patient Records - Microsoft Internet Explorer Links >> AHCCCS Health Information Exchange Viewer






View Record Detail






View Record Detail






Navigation: Return to Record List






Multiple Pop-ups Open





HIE Messages



Clinical System Gateway &

RLS Gateway
Publish

Update

RLS Gateway

&

Secure Web Viewer Gateway

Patient Lookup/Search

Search Results

Data Listing Request

Data Listing Response

Secure Web Viewer Gateway

&

Clinical System Gateway

Data Request Data Response



HIeHR Utility Web Services

 Clinical System Gateway and RLS Gateway

- Publish
 - Clinical System informs RLS of data record
 - meta-data, identifies patient, indicates type of record
- Update
 - Indicates Changes to previously published records



HIeHR Utility Web Services

- RLS Gateway and Secure Web Viewer Gateway
 - Patient Lookup/Search
 - Secure Web Viewer Gateway submits search criteria to RLS Gateway
 - Patient Search Results
 - RLS Gateway returns matching patients to Secure Web Viewer Gateway
 - Data List Request
 - Secure Web Viewer Gateway submits request for patient data record listing to RLS Gateway
 - Data Listing Response
 - RLS Gateway returns patient data record listing to Secure Web Viewer Gateway




HIeHR Utility Web Services

 Secure Web Viewer Gateway and Clinical System Gateway


- Data Request
 - Secure Web Viewer Gateway submits request for specific data record to Clinical System Gateway
- Data Response
 - Clinical System Gateway returns specific data record to Secure Web Viewer Gateway



HIE Message Structure



HL 7 3.0
Publish
Update
Lookup
Data Request
Patient Search
Record Request



CCD
Discharge Summary
Lab Results
Medication History




AHCCCS Strategy for Using CCD Levels of Encoding

- Level of data encoding may vary among **AHCCCS HIE participants**
- Custom markup can augment CCD markup
- AHCCCS will provide
 - data partner implementation guide
 - business rules
 - data elements
 - xml



CCD Structure

- XML document
 - header
 - body
- Header identifies patient, provider, document type, originating organization
- Body has a
 - human-readable section
 - Human-readable part is attested and contains complete content
 - Either XML or any MIME encoded BLOB e.g. *.doc, *.pdf, *.tif
 - encoded section
 - if present, can be safely ignored by recipients unable to process it



SHO

Header

- Document Type
- Sender
- Receiver
- Patient

Body

T E X

Section(s)

- Admission Details
- Primary/Secondary Diagnosis
- Observations
- Medications
- •Follow-up

0

Entries

- Admission Details
- Primary/Secondary Diagnosis
- Observations
- Medications



Sample components of CCD

```
<ClinicalDocument>
                                                 Header
      <structuredBody>
            <section>
                  <text>...</text>
                                         Narrative Block
                  <observation>...
                  <substanceAdministration>
                          <supply>...</supply>
                  </substanceAdministration>
            </section>
            <section>...</section>
      </structuredBody>
</ClinicalDocument>
```




Example of a fully encoded document section

```
<section>
  <title>Vital Signs</title>
  <text>Temperature is 36.9 C</text>
  <entry>
     <observation classCode="OBS" moodCode="EVN">
        <code code="386725007" codeSystem="2.16.840.1.113883.6.96"</pre>
           codeSystemName="SNOMED CT" displayName="Body temperature"/>
        <effectiveTime value="200004071430"/>
        <value xsi:type="PQ" value="36.9" unit="Cel"/>
     </observation>
  </entry>
</section>
```




Clinical System and Gateway Communications






Clinical System and Gateway Communications






Clinical System and Gateway Communications





RLS Architecture






'Defense in Depth' Security

Application Based Security User Authentication Virus Scanning

Server Configuration Management
User Access Management
Database Encryption

Layer 3





Our first care is your healthcare